

## DPM microdoser



Moretto S p A offers high-performance and precise 'DPM' microdoser. This microdoser is a microbatch pulsed blender that interprets the dosage. It consists of a rotopulse device and double tower that are the exclusive equipment of Moretto. These devices allow dosing a very small quantity of the product with a microbatch of 0.1 gm, precision of 0.001 per cent and reaction time of 25 ms. This microdoser is unique because it exploits a static system to determine weights. It first doses master and then determines higher percentages for the other components through a series of 'micropulsations' based on the small product quantity. It is a pulsed dosing unit with direct setting of the percentages. It does not need a mixer as it simultaneously doses all components; moreover, all its hoppers can be removed. The microdoser can store up to 100 recipes and does not accumulate material, thus ensuring the easy change of production. The reduced dimensions allow the direct installation on machine board. It is available in four basic models: 'DPM 0', 'DPM 1', 'DPM 2' and 'DPM 4' and allows productions up to four components.

### **Moretto S p A**

Via dell'Artigianato, 3  
35010, Massanzago/PD Italia.

Tel: +39-049-939 6711

Fax: +39-049-939-6710

Email: [marketing@moretto.com](mailto:marketing@moretto.com)